

Product Datasheet

p57 Kip2 (Phospho-T310) antibody (orb669866)

Description	Rabbit polyclonal antibody to p57 Kip2 (Phospho-T310)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic phosphopeptide derived from human p57 Kip2 around the phosphorylation site of Threonine 310.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 50 kDa
Uniprot ID	P49918
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.