


Product Datasheet

KDM5C antibody (orb66993)

Description	Rabbit polyclonal to KDM5C. Lysine-specific demethylase 5C (KDM5C) is a histone demethylase that spec...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human KDM5C
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	176
Uniprot ID	P41229
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...