


Product Datasheet

KDM3A antibody (orb67000)

Description	Rabbit polyclonal to KDM3A. Lysine-specific histone demethylase 3A (KDM3A) is a histone demethylase t...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human KDM3A
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	147
Uniprot ID	Q9Y4C1
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...