

Product Datasheet

PERFORIN antibody (orb67011)

Description	Rabbit polyclonal to PERFORIN.PERFORIN (PRF1) is a cytolytic protein found in the granules of CD8 T-...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of mouse PERFORIN
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	62
Uniprot ID	P31240
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.

