


Product Datasheet

EOMES antibody (orb67012)

Description	Rabbit polyclonal to EOMES.Eomesodermin homolog (EOMES) functions as a transcriptional activator pla...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of mouse EOMES
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	75
Uniprot ID	O54839
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...