



Product Datasheet

PIAS1 antibody (orb67014)

Description	Rabbit polyclonal to PIAS1.E3 SUMO-protein ligase PIAS1 is a member of the mammalian PIAS (protein i...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, IP, WB
Immunogen	C-terminal region of human PIAS1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	antiserum
MW	72
Uniprot ID	O88907
Dilution Range	WB:1:1000-1:5000,IP:1:50-1:100,IF:1:100-1:200,IHC:1:100-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...



Western blot analysis of total cell extr...