


Product Datasheet

CARM1 antibody (orb67015)

Description	Rabbit polyclonal to CARM1. Coactivator-associated arginine methyltransferase 1 (CARM1) is also known...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human CARM1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	66
Uniprot ID	Q86X55
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...