

Product Datasheet

GPR180 antibody (orb670162)

Description	Rabbit polyclonal antibody to GPR180
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide, corresponding to Human POTE _G .
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.4.
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -25°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 58,54 kDa
Uniprot ID	Q6S545, Q6S5H5, A6NI47
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.