

Product Datasheet

GPR196 antibody (orb670178)

Description	Rabbit polyclonal antibody to GPR196
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide, corresponding to Human RGPD1/2/3/4/5/8.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.4.
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -25°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 198 kDa
Uniprot ID	P0DJD0 , Q7Z3J3 , Q99666 , A6NKT7 , P0DJD1 , O14715
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.