


Product Datasheet

WNT2 antibody (orb67019)

Description	Rabbit polyclonal to WNT2.WNT2 belongs to the WNT family, which is a ligand for members of the frizz...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of mouse WNT2
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	40
Uniprot ID	P21552
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...