


Product Datasheet

GRP78 antibody (orb67022)

Description	Rabbit polyclonal to GRP78.78 kDa glucose-regulated protein (GRP78) belongs to the heat shock protei...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of mouse GRP78
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	72
Uniprot ID	P20029
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...