

Product Datasheet

PGP9.5 antibody (orb67029)

Description	Rabbit polyclonal antibody to PGP9.5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human UCHL1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	25
Uniprot ID	P09936
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.

