


Product Datasheet

MMS21 antibody pair (orb67033)

Description	Rabbit polyclonal to MMS21. 50 μg MMS21 (N-term)&50 μg MMS21 (C-term). vPair antibodies repres...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	-
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	28
Uniprot ID	Q96MF7
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...