


Product Datasheet

DNMT2 antibody (orb67038)

Description	Rabbit polyclonal to DNMT2. DNA (cytosine-5)-methyltransferase-like protein 2 (DNMT2) is also known as...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human DNMT2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	45
Uniprot ID	O14717
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...