


Product Datasheet

FOXP1 antibody (orb67039)

Description	Rabbit polyclonal to FOXP1. Forkhead box protein P1 (FOXP1) is a transcriptional repressor that plays...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human FOXP1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	75
Uniprot ID	Q9H334
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...