


Product Datasheet

FZD8 antibody (orb67041)

Description	Rabbit polyclonal to FZD8. Frizzled-8 (FZD8) is a receptor for Wnt proteins, which belongs to the G-p...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human FZD8
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	73
Uniprot ID	Q9H461
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...