


Product Datasheet

GAPDH antibody (orb670411)

Description	GAPDH antibody
Species/Host	Rabbit
Reactivity	Gallus, Human, Monkey, Mouse, Other, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant Protein
Target	GAPDH
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04406
Dilution Range	WB: 1: 5000IHC-P: 1: 200
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T (Lane 1), ...