

Product Datasheet

ITGA4 antibody (orb67043)

Description	Rabbit polyclonal to ITGA4. Integrin alpha-4 (ITGA4) is a single-pass type I membrane protein that be...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human ITGA4
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	115
Uniprot ID	P13612
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.

