


Product Datasheet

Acetyl-Histone H4 (Lys16) antibody (orb670456)

Description	Acetyl-Histone H4 (Lys16) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	HIST1H4A
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62805
Dilution Range	WB: 1: 1000-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Hela cells using...