


Product Datasheet

RSF1 antibody (orb67050)

Description	Rabbit polyclonal to RSF1. Remodeling and spacing factor 1 (SF1) is a nuclear protein required for ass...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human RSF1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	164
Uniprot ID	Q96T23
Dilution Range	WB: 1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...