


Product Datasheet

HDAC9 antibody (orb67059)

Description	Rabbit polyclonal to HDAC9. Histone deacetylase 9 (HDAC9) belongs to the histone deacetylase family, ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human HDAC9
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	111
Uniprot ID	Q9UKV0
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...