


Product Datasheet

SMARCA5 antibody (orb67068)

Description	Rabbit polyclonal to SMARCA5.SWI/SNF-related matrix-associated actin-dependent regulator of chromati...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human SMARCA5
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	122
Uniprot ID	O60264
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...