


Product Datasheet

LAT antibody (orb67078)

Description	Rabbit polyclonal to LAT.Linker for activation of T-cells family member 1 (LAT) is a single-pass typ...
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of mouse LAT
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	26
Uniprot ID	O54957
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...