


Product Datasheet

HSP90AA1 antibody (orb67079)

Description	Rabbit polyclonal to HSP90AA1. Heat shock protein HSP 90 (HSP90) proteins are highly conserved molecu...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human HSP90AA1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	85
Uniprot ID	P07900
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...