


Product Datasheet

CD20 antibody (orb67092)

Description	Rabbit polyclonal to CD20.B-lymphocyte antigen CD20 is a member of the membrane-spanning 4A gene fam...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human CD20
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	33
Uniprot ID	P11836
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...