


Product Datasheet

IL2RB antibody (orb67097)

Description	Rabbit polyclonal to IL2RB. Interleukin-2 receptor subunit beta (IL2RB) is a type I membrane protein ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of mouse IL2RB
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	61
Uniprot ID	P16297
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...