

Product Datasheet

BACH2 antibody (orb67100)

Description	Rabbit polyclonal to BACH2. Transcription regulator protein BACH2 belongs to the bZIP family and the ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human BACH2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	93
Uniprot ID	Q9BYV9
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.

