


Product Datasheet

FOS antibody (orb67101)

Description	Rabbit polyclonal to FOS. The Fos gene family consists of 4 members: FOS, FOSB, FOSL1, and FOSL2. The...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of mouse FOS
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	41
Uniprot ID	P01101
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...