


Product Datasheet

HPRT1 antibody (orb67108)

Description	Rabbit polyclonal to HPRT1. Hypoxanthine-guanine phosphoribosyltransferase (HPRT1) belongs to the pur...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human HPRT1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	25
Uniprot ID	P00492
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...