


Product Datasheet

NOTCH2 antibody (orb67110)

Description	Rabbit polyclonal antibody to NOTCH2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human NOTCH2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	265
Uniprot ID	Q04721
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extracellular NOTCH2.