


Product Datasheet

HOXA9 antibody (orb67115)

Description	Rabbit polyclonal to HOXA9. In vertebrates, the genes encoding the class of transcription factors cal...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human HOXA9
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	30
Uniprot ID	P31269
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...