


Product Datasheet

ADCY1 antibody (orb67118)

Description	Rabbit polyclonal to ADCY1. Adenylate cyclase type 1 (ADCY1) is a multi-pass membrane protein that be...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	N-terminal region of human ADCY1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	123
Uniprot ID	Q08828
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot of total cell extracts from...