


Product Datasheet

SYT5 antibody (orb67119)

Description	Rabbit polyclonal to SYT5. Synaptotagmins, such as SYT5, are a family of type III membrane proteins c...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human SYT5
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	43
Uniprot ID	O00445
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...