


Product Datasheet

AXIN2 antibody (orb67132)

Description	Rabbit polyclonal to AXIN2. Axin-2 (AXIN2) may act as an inhibitor of the Wnt signaling pathway by do...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human AXIN2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	94
Uniprot ID	Q9Y2T1
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...