

Product Datasheet

SOX9 antibody (orb67136)

Description	Rabbit polyclonal to SOX9. Transcription factor SOX-9 is a nucleic protein that plays an important ro...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human SOX9
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	56
Uniprot ID	P48436
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.

