

Product Datasheet

OR7C1 antibody (orb671378)

Description	OR7C1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 7C1.
Target	OR7C1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	076099
Dilution Range	WB: 1:500-1:2000, IF: 1:200-1:1000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.

