


## Product Datasheet

### DKK1 antibody (orb67138)

|                            |                                                                                                          |
|----------------------------|----------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal to DKK1. Dickkopf-related protein 1 (DKK1) belongs to the dickkopf family. DKK1 is ... |
| <b>Species/Host</b>        | Rabbit                                                                                                   |
| <b>Reactivity</b>          | Human                                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                                             |
| <b>Tested Applications</b> | WB                                                                                                       |
| <b>Immunogen</b>           | N-terminal region of human DKK1                                                                          |
| <b>Preservatives</b>       | 0.01% NaN <sub>3</sub>                                                                                   |
| <b>Storage</b>             | -20°C                                                                                                    |
| <b>Note</b>                | For research use only                                                                                    |
| <b>Clonality</b>           | Polyclonal                                                                                               |
| <b>Purity</b>              | Affinity purification                                                                                    |
| <b>MW</b>                  | 27                                                                                                       |
| <b>Uniprot ID</b>          | <a href="https://www.uniprot.org/entry/O94907">O94907</a>                                                |
| <b>Dilution Range</b>      | WB:1:500-1:1000                                                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                          |



Western blot analysis of total cell extr...