


Product Datasheet

WHSC1 antibody (orb67150)

Description	Rabbit polyclonal to WHSC1. Wolf-Hirschhorn syndrome candidate 1 protein (WHSC1) is a probable histon...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human WHSC1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	152
Uniprot ID	O96028
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...