


Product Datasheet

APOBEC3 gamma antibody (orb67157)

Description	Rabbit polyclonal to APOBEC3 gamma.DNA dC-dU-editing enzyme APOBEC-3G (APOBEC3G) belongs to the cyti...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human APOBEC3G
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	46
Uniprot ID	Q9HC16
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot of total cell extracts from...