


Product Datasheet

ATP11A antibody (orb67158)

Description	Rabbit polyclonal to ATP11A. Probable phospholipid-transporting ATPase 1H (ATP11A) is a multi-membran...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human ATP11A
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	130
Uniprot ID	P98196
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...