


Product Datasheet

CBL antibody (orb67165)

Description	Rabbit polyclonal to CBL.E3 ubiquitin-protein ligase CBL is an adapter protein that functions as a n...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human CBL
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	100
Uniprot ID	P22681
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...