


Product Datasheet

CARD10 antibody (orb67173)

Description	Rabbit polyclonal to CARD10. The caspase recruitment domain (CARD) is a protein module that consists of...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human CARD10
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	116
Uniprot ID	Q9BWT7
Dilution Range	WB: 1:1000-1:3000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...