


Product Datasheet

DCAMKL1 antibody (orb67177)

Description	Rabbit polyclonal to DCAMKL1. Serine/threonine-protein kinase DCLK1 (DCAMKL1) belongs to the protein ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human DCAMKL1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	82
Uniprot ID	O15075
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...