


Product Datasheet

E4F1 antibody (orb67179)

Description	Rabbit polyclonal to E4F1. Transcription factor E4F1 may function as a transcriptional repressor. It ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human E4F1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	83
Uniprot ID	Q66K89
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...