


## Product Datasheet

### GTF2IRD1 antibody (orb671819)

|                            |                                                                              |
|----------------------------|------------------------------------------------------------------------------|
| <b>Description</b>         | GTF2IRD1 antibody                                                            |
| <b>Species/Host</b>        | Rabbit                                                                       |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                 |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                               |
| <b>Immunogen</b>           | Synthesized peptide derived from the N-terminal region of Human WBSCR11.     |
| <b>Target</b>              | GTF2IRD1                                                                     |
| <b>Preservatives</b>       | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide. |
| <b>Form/Appearance</b>     | Liquid                                                                       |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                |
| <b>Note</b>                | For research use only                                                        |
| <b>Isotype</b>             | IgG                                                                          |
| <b>Clonality</b>           | Polyclonal                                                                   |
| <b>Uniprot ID</b>          | <a href="#">Q9UHL9</a>                                                       |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IHC-P: 1:100-1:300, ELISA: 1:20000                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                              |



Western blot analysis of 293 cells using...