

Product Datasheet

FBW7 antibody (orb67185)

Description	Rabbit polyclonal to FBW7.F-box/WD repeat-containing protein 7 (FBXW7) is a substrate recognition co...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human DDB1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	80
Uniprot ID	Q969H0
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.

