


Product Datasheet

IL13 Receptor alpha 2 antibody (orb67188)

Description	Rabbit Polyclonal antibody to IL13 Receptor alpha 2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human IL13RA2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	44
Uniprot ID	Q14627
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...