

Product Datasheet

CCR5 antibody (orb67199)

Description	Rabbit polyclonal to CCR5.C-C chemokine receptor type 5 (CCR5) is a multi-pass membrane protein that...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human CCR5
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	41
Uniprot ID	P51681
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.

