


Product Datasheet

CD80 antibody (orb67204)

Description	Rabbit polyclonal to CD80.T-lymphocyte activation antigen CD80 is a single-pass type I membrane prot...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human CD80
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	33
Uniprot ID	P33681
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...