


Product Datasheet

CX3CR1 antibody (orb67205)

Description	Rabbit polyclonal to CX3CR1. CX3C chemokine receptor 1 (CX3CR1) is a multi-pass membrane protein that...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human CX3CR1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	40
Uniprot ID	P49238
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...